

# Dipslides

Simple and economic to use, Palintest Dipslides are ideal for cooling water circuits, institutional hygiene and applications where indication of general microbiological contamination is required.

## Palintest Dipslides include the following:

- **Suitable for liquid and surface testing** – either immerse the Dipslide within the sample or place in contact with the surface under test.

- **Available in two formats** – supplied in packs of 10, Dipslides are suitable for Total Colony (TTC) only or Total Colony plus Yeasts and Mould (TTC/M) for full compliance checking.

- **Simple result determination** – compare the clearly visible colonies produced with the supplied standard chart.

Palintest Dipslides are suitable for use within the range of incubators on page 62.



Ordering Information	Part Code
Dipslides, TTC (10 pack)	PT711
Dipslides, TTC/M (10 pack)	PT710